


GOLDEN 3 24, 26, 28 - diameters of lines COUSIN [mm]

MIDDLE OF THE CANNOPY



- 0.6 = dyneema LIROS DC 060
- 0.8 = dyneema LIROS DC 120
-
- 0.95 =dyn. COUSIN 0.95 mm
- 1.1 = dyn. COUSIN 1.1 mm
- 1.3 = dyn. COUSIN 1.3 mm
-
- 1.9 = Technora COUSIN 1.9
- 1.8 = Technora COUSIN 1.8
- 1.5 = Technora COUSIN 1.5
- 1.2 = Technora COUSIN 1.2